HPLC 用充填カラム Develosil ANIDIUS

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PACKED COLUMN Develosil ANIDIUS INSTRUCTION MANUAL

Introduction

Thank you for purchasing our Develosil ANIDIUS. We recommend you to read this manual carefully to ensure maximal efficiency and long life of the column.

Characteristic of Develosil ANIDIUS

Develosil ANIDIUS was designed for Hydrophilic Interaction Liquid Chromatography (HILIC). Develosil ANIDIUS is most suitable for the analysis of the high polarity compounds which can't maintain it by the reverse phase chromatography.

Characteristic of Develosil ANIDIUS

Silica base / High purity Silica gel, Particle size / 5 μ m, Surface area / 450m²/g

Pore diameter / 14nm, pore capacity / 1.15ml/g, Bonded phase / -NH₂, -OH

Care and Use

- In general, acetonitrile/water is used as mobile phase.
- When mobile phase has high acetonitrile content, note dissolubility of the salt. The dissolubility of phosphate buffers used widely in reserved phase chromatography is low in acetonitrile. Therefore use of phosphate buffers is not recommended.
- Develosil ANIDIUS is designed for use in the pH range of 2-8.
- Whenever possible, samples should be dissolved in a solvent that has the same composition as the mobile phase.
- Keep column pressure under 20MPa.
- Use acetonitrile content in the mobile phase within the range of 40-95%.
- Please read 「**DEVELOASIL PACKED COLUMN INSTRUCTION MANUAL**」 before the use.

Develosil Columns are manufactured, inspected, packaged and shipped under strict standards of quality control. Should you find any defect in performance, please contact us. However, that we can't guarantee the lifetime of columns, nor can we accept any claim when their performance has deteriorated due to noncompliance with the above operating instructions.

NOMURA CHEMICAL CO., LTD.

15 Hinode-cho, Seto, Aichi 489-0004 Japan

TEL +81-561-48-1853 FAX +81-561-48-1434

E-mail info@develosil.net